



Department of Physics
Rishi Bankim Chandra Evening College
 Naihati, North 24 Parganas, PIN – 743165, West Bengal, India
 Phone: (033) 2581-1281 ♦ Email: physics@rbcec.in

Determination of Modulus of Rigidity (n) of the Material of a Wire by Dynamical Method

- Theory:**




Figure: Instrument for determining rigidity modulus

The time-period T with which the bob of a torsion pendulum oscillates, with its suspension wire as axis, is given by,

$$T = 2\pi \sqrt{I/c}$$

$$\text{Or, } c = \frac{4\pi^2 I}{T^2} \dots (1)$$

where, I is the moment of inertia of the suspended cylinder about its own axis and is given by,

$$I = \frac{1}{2} \times M \times R^2 \dots (2)$$

where, M = mass of the bob

R = radius of the bob

and c represents the couple exerted by the suspension wire of length l and radius r for one radian twist at its free end and is given by,

$$c = \frac{\eta \pi r^4}{2l} \quad \dots (3)$$

where, n = rigidity modulus of the material of the wire.

From equations (1) and (3), we can write,

$$\frac{4\pi^2 I}{T^2} = \frac{\eta \pi r^4}{2l}$$

$$\text{Or, } \eta = \frac{8\pi l}{T^2 r^4} \times I \quad \dots (4)$$

At first calculate I from equation (2). Then by measuring r , l and T calculate η from equation (4). SI unit of rigidity modulus will be in N/m^2 .

(next page)

• **Experimental Data:**

(A) Determination of the moment of inertia (I) of the suspended cylinder:

(i) Mass of the cylinder = $M = \dots\dots\dots$ kg (given)

(ii) Determination of Vernier constant:

1 small division of main scale = 0.1 cm

Now, 10 vernier divisions = 9 small main scale divisions

$$\text{or, } 1 \text{ v.d.} = \frac{9}{10} \text{ s.m.d}$$

So, Vernier constant = 1 small main scale division – 1 vernier division

$$= \left(1 - \frac{9}{10}\right) \times 0.1 = 0.01 \text{ cm}$$

(iii) Determination of mechanical error (e) of Slide Callipers:

$$e = \pm y \times v \quad \text{cm} \quad [\text{where } y \text{ is the initial reading of erroneous slide calipers}]$$

Table – I: Calculation of Radius (R) of the Cylinder

No. of Observations	Measuring Quantity (Diameter)	Main Scale Reading (In cm)	Vernier Reading	Total Reading = $S + N \times v$	Mean D' (In cm)	Corrected Diameter $D = D' - e$ (In cm)	Radius (R) (In m)
1	Diameter at one point						
	Diameter at another point perpendicular to previous						
2							
3							
4							
5							

So from equation (2),

$$I = \frac{1}{2} \times M \times R^2 = \dots\dots\dots \text{kg-m}^2$$

(B) Measurement of the Length (l) of the suspension wire:

$$l = \frac{\dots + \dots + \dots}{3} \text{ m} = \dots\dots\dots \text{ m}$$

(C) Measurement of the radius (r) of the suspension wire:

Measuring Least Count of Screw Gauge:

Smallest division of least count = S = mm

Pitch of the screw = p = mm

Number of divisions on the circular scale = N =.....

Least count of instrument = L.C. = $\frac{P}{N} = \dots\dots \text{ mm}$

mechanical error (e) of Screw Gauge:

$$e = \pm y \times v \quad \text{mm} \quad [\text{where } y \text{ is the initial reading of erroneous screw gauge}]$$

Table – II: Calculation of Radius (r) of the Suspension Wire

Serial No.	Measuring Quantity (Diameter)	Linear Scale Reading (LSR) (In mm)	Circular Scale Reading (CSR)	Diameter (in mm) Total = LSR+ (C.S.R X L.C.)	Mean D' (In mm)	Corrected Diameter D = D' - e (In mm)	Radius (r) (In m)
1	Diameter at one point						
	Diameter at another point perpendicular to previous						
2							
3							
4							
5							

(D) Determination of the time-period of torsional oscillation (T) of the cylinder:

Table – III

No. of obs.	Time for 30 oscillations (in s)	Mean time <i>t</i> (in s)	Time-period $T = \frac{t}{30}$ (in s)
1.
2.		
3.		
4.		
5.		

• Calculations:

From experiment, we have,

$$I = \dots\dots\dots \text{kg-m}^2$$

$$l = \dots\dots\dots \text{m}$$

$$r = \dots\dots\dots \text{m}$$

$$T = \dots\dots\dots \text{s}$$

Then, $\eta = \frac{8\pi l}{T^2 r^4} \times I = \dots\dots\dots \text{N/m}^2$